

Product Information	
Catalog No.	PA10146
Product name	Anti-AGR2 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Pongo abelii
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	AGR2, AG-2, AG2, Anterior gradient protein 2 homolog, hAG-2, HPC8, Secreted cement gland protein XAG-2 homolog
Immunogen	E. coli - derived recombinant Human AGR2 (Arg21-Leu175).
Accession	O95994
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.