

Product Information	
Catalog No.	PA10153
Product name	Anti-ALDH1A1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat, Bos taurus, Equus caballus, Macaca fascicularis, Macroscelides proboscideus, Mesocricetus auratus, Oryctolagus cuniculus, Ovis aries
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Aldehyde dehydrogenase, cytosolic, Retinal dehydrogenase 1, ALDH-E1, RalDH1, Aldehyde dehydrogenase family 1 member A1, ALDH1A1, ALDH1, RALDH 1, ALDC, ALHDII, PUMB1
Immunogen	E. coli - derived recombinant Human ALDH1A1 (Met1-Ser501).
Accession	P00352
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.