

Product Information	
Catalog No.	PA10171
Product name	Anti-CA1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Gorilla gorilla gorilla, Macaca mulatta, Macaca nemestrina, Pan troglodytes
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	CA-I, Carbonate dehydratase I, Carbonic anhydrase B, CA1, CAB, Carbonic anhydrase 1, Carbonic anhydrase I
Immunogen	E. coli - derived recombinant Human CA1 (Ala2-Phe261).
Accession	P00915
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.