

Product Information	
Catalog No.	PA10181
Product name	Anti-MYC Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Cat, Rat, Bos taurus, Callithrix jacchus, Dog, Feline leukemia virus FTT, Feline leukemia virus, Galeopterus variegatus, Hylobates lar, Macaca mulatta, Marmota monax, Ovis aries, Pan paniscus, Pan troglodytes, Pongo pygmaeus, Pteropus hypomelanus, Pig, Sylvilagus floridanus, Tadarida brasiliensis, Tupaia glis
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Proto-oncogene c-Myc, BHLHE39, MYC, Myc proto-oncogene protein, bHLHe39, Transcription factor p64, Class E basic helix-loop-helix protein 39
Immunogen	E. coli - derived recombinant Human MYC (Glu351-Ser437).
Accession	P01106
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.