

Product Information	
Catalog No.	PA10190
Product name	Anti-AVP Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Balaenoptera physalus, Bos taurus, Cavia porcellus, Equus caballus, Human, Loxodonta africana, Ovis aries, Pig
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Neurophysin-II, Arginine-vasopressin, VP, ARVP, AVP-NP11, AVP, Vasopressin-neurophysin 2-copeptin
Immunogen	E. coli - derived recombinant Human AVP (Cys20-Tyr164).
Accession	P01185
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.