

Product Information	
Catalog No.	PA10191
Product name	Anti-TSHB Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Cat, Rat, Bos taurus, Dog, Equus caballus, Lama glama, Phodopus sungorus, Pig
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Thyrotropin subunit beta, TSHB, Thyrotropin beta chain, Thyroid-stimulating hormone subunit beta, TSH-beta, TSH-B, Thyrotropin alfa
Immunogen	E. coli - derived recombinant Human TSHB (Phe21-Tyr132).
Accession	P01222
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.