

Product Information	
Catalog No.	PA10238
Product name	Anti-ESR1/ER-alpha Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Cat, Rat, Anolis carolinensis, Aspidoscelis uniparens, Bos taurus, Equus caballus, Chicken, Macaca mulatta, Mesocricetus auratus, Pig, Taeniopygia guttata
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Estrogen receptor, ER-alpha, ESR, ER, ESR1, Nuclear receptor subfamily 3 group A member 1, NR3A1, Estradiol receptor
Immunogen	E. coli - derived recombinant Human ESR1/ER-alpha (Met297-Ser554).
Accession	P03372
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.