

Product Information	
Catalog No.	PA10272
Product name	Anti-CYP17A1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Macaca fascicularis, Macaca mulatta, Pan troglodytes, Papio cynocephalus, Papio hamadryas ursinus
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Steroid 17-alpha-hydroxylase/17, 20 lyase, CYPXVII, Cytochrome P450 17A1, Cytochrome P450-C17, Steroid 17-alpha-monooxygenase, 17-alpha-hydroxyprogesterone aldolase, CYP17A1, S17AH, Cytochrome P450c17, CYP17
Immunogen	E. coli - derived recombinant Human CYP17A1 (Cys183-Thr508).
Accession	P05093
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.