

Product Information	
Catalog No.	PA10275
Product name	Anti-IL4 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Cercopithecus atys, Macaca fascicularis, Macaca mulatta, Pan troglodytes, Papio anubis
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Binetrakin, B-cell stimulatory factor 1, Pitratkinra, Lymphocyte stimulatory factor 1, IL-4, BSF-1, IL4, Interleukin-4
Immunogen	E. coli - derived recombinant Human IL4 (Gly24-Ser153).
Accession	P05112
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.