

Product Information	
Catalog No.	PA10311
Product name	Anti-LDHB Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Mouse, Alligator mississippiensis, Anas platyrhynchos, Bos taurus, Caiman crocodilus apaporiensis, Chicken, Human, Macaca fascicularis, Monodelphis domestica, Oryctolagus cuniculus, Pan troglodytes, Pelodiscus sinensis japonicus, Rat, Sceloporus woodi, Pig, Trachemys scripta elegans
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	LDHB, LDH-B, L-lactate dehydrogenase B chain, LDH heart subunit, LDH-H, Renal carcinoma antigen NY-REN-46
Immunogen	E. coli - derived recombinant Human LDHB (Met1-Leu334).
Accession	P07195
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.