


Product Information	
Catalog No.	PA10312
Product name	Anti-NEFL Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat, Bos taurus, Pongo abelii, Pig
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	NEFL, Neurofilament light polypeptide, 68 kDa neurofilament protein, NF68, Neurofilament triplet L protein, NF-L, NFL
Immunogen	E. coli - derived recombinant Human NEFL (Met1-Asn352).
Accession	P07196
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-NEFL polyclonal antibody at 1ug/mL on various samples.

- Lane 1: Mouse brain lysate
- Lane 2: Mouse cerebellum lysate
- Lane 3: Rat brain lysate
- Lane 4: Rat cerebellum lysate