


Product Information	
Catalog No.	PA10350
Product name	Anti-CD41/ITGA2B Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	CD41, ITGAB, GPIIb, ITGA2B, Platelet membrane glycoprotein IIb, GP2B, Integrin alpha-IIb, GPalpha IIb
Immunogen	E. coli - derived recombinant Human CD41/ITGA2B (Ala678-Leu894).
Accession	P08514
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-ITGA2B polyclonal antibody at 1ug/mL on various samples.
Lane 1: Mouse spleen lysate
Lane 2: Mouse lung lysate
Lane 3: Rat spleen lysate
Lane 4: Rat lung lysate