

Product Information	
Catalog No.	PA10451
Product name	Anti-MYH1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Dog, Rat, Bos taurus, Equus caballus, Chicken, Oryctolagus cuniculus, Pig
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Myosin heavy chain 2x, Myosin heavy chain 1, Myosin-1, MyHC-IIx/d, MyHC-2x, Myosin heavy chain, skeletal muscle, adult 1, Myosin heavy chain IIx/d, MYH1
Immunogen	E. coli - derived recombinant Human MYH1 (Pro842-Lys967).
Accession	P12882
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.