


Product Information	
Catalog No.	PA10452
Product name	Anti-MYH7 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Dog, Rat, Bos taurus, Equus caballus, Mesocricetus auratus, Pig
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	MYH7, Myosin heavy chain slow isoform, MyHC-beta, MyHC-slow, MYHCB, Myosin heavy chain 7, Myosin heavy chain, cardiac muscle beta isoform, Myosin-7
Immunogen	E. coli - derived recombinant Human MYH7 (Arg1268-Leu1516).
Accession	P12883
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot analysis was performed using anti-MYH7 polyclonal antibody at 1ug/mL on various samples.

Lane 1: Mouse heart lysate

Lane 2: Mouse skeletal muscle lysate

Lane 3: Rat skeletal muscle lysate

Western blot