

Product Information	
Catalog No.	PA10471
Product name	Anti-MIF Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Bos taurus, Human, Macaca fascicularis, Macaca mulatta, Meriones unguiculatus, Mouse, Ovis aries, Rat, Pig
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	MIF, L-dopachrome isomerase, L-dopachrome tautomerase, GIF, Macrophage migration inhibitory factor, MMIF, Phenylpyruvate tautomerase, Glycosylation-inhibiting factor, GLIF
Immunogen	E. coli - derived recombinant Human MIF (Met1-Ala115).
Accession	P14174
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.