


Product Information	
Catalog No.	PA10502
Product name	Anti-NQO1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Pongo abelii, Rat
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	NAD(P)H:quinone oxidoreductase 1, Azoreductase, Phylloquinone reductase, NQO1, DT-diaphorase, NAD(P)H dehydrogenase [quinone] 1, QR1, DIA4, NMOR1, Quinone reductase 1, Menadione reductase, DTD
Immunogen	E. coli - derived recombinant Human NQO1 (Met1-Lys274).
Accession	P15559
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-NQO1 polyclonal antibody at 1ug/mL on various samples.

Lane 1: Lncap cell lysate

Lane 2: LO2 cell lysate

Lane 3: A549 cell lysate