

Product Information	
Catalog No.	PA10537
Product name	Anti-PSMC3 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat, Xenopus laevis
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	26S proteasome AAA-ATPase subunit RPT5, Tat-binding protein 1, Proteasome 26S subunit ATPase 3, Proteasome subunit P50, 26S proteasome regulatory subunit 6A, TBP-1, PSMC3, TBP1
Immunogen	E. coli - derived recombinant Human PSMC3 (Asn2-Ala439).
Accession	P17980
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.