


Product Information	
Catalog No.	PA10592
Product name	Anti-BMP6 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	VG-1-R, VG-1-related protein, Bone morphogenetic protein 6, VGR-1, BMP-6, VGR, BMP6
Immunogen	E. coli - derived recombinant Human BMP6 (Ser375-His513).
Accession	P22004
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot analysis was performed using anti-BMP6 polyclonal antibody at 1ug/mL on various samples.
Lane 1: Mouse heart lysate

Western blot