

Product Information	
Catalog No.	PA10598
Product name	Anti-IGFBP4 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat, Bos taurus, Ovis aries, Pig
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	IGF-binding protein 4, IBP-4, Insulin-like growth factor-binding protein 4, IGFBP-4, IGFBP4, IBP4
Immunogen	E. coli - derived recombinant Human IGFBP4 (Asp22-Glu258).
Accession	P22692
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.