

Product Information	
Catalog No.	PA10679
Product name	Anti-GDI1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Bos taurus, Dog, Human, Macaca fascicularis, Mouse, Pan troglodytes, Pongo abelii, Pongo pygmaeus, Rat, Pig
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Oligophrenin-2, Rab GDP dissociation inhibitor alpha, XAP4, Guanosine diphosphate dissociation inhibitor 1, GDIL, GDI1, RABGDIA, Rab GDI alpha, Protein XAP-4, OPHN2, GDI-1
Immunogen	E. coli - derived recombinant Human GDI1 (Met1-Val274).
Accession	P31150
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.