

Product Information	
Catalog No.	PA10702
Product name	Anti-MCM5 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Bos taurus, Caenorhabditis elegans, Drosophila melanogaster, Human, Mouse, Xenopus tropicalis, Xenopus laevis
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	CDC46, DNA replication licensing factor MCM5, CDC46 homolog, P1-CDC46, MCM5
Immunogen	E. coli - derived recombinant Human MCM5 (Val331-Ile537).
Accession	P33992
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.