

Product Information	
Catalog No.	PA10719
Product name	Anti-PTGS2 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat, Bos taurus, Cavia porcellus, Equus caballus, Chicken, Neovison vison, Oryctolagus cuniculus, Ovis aries
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	PGH synthase 2, Prostaglandin G/H synthase 2, Cyclooxygenase-2, PTGS2, COX2, PGHS-2, PHS II, Prostaglandin-endoperoxide synthase 2, COX-2, Prostaglandin H2 synthase 2
Immunogen	E. coli - derived recombinant Human PTGS2 (Phe187-Asn425).
Accession	P35354
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.