

Product Information	
Catalog No.	PA10723
Product name	Anti-MYH10 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Dog, Rat, Bos taurus, Chicken, Pig
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	NMMHC II-b, MYH10, Non-muscle myosin heavy chain IIb, NMMHC-IIb, Cellular myosin heavy chain, type B, Myosin heavy chain 10, Myosin heavy chain, non-muscle IIb, Myosin-10, Non-muscle myosin heavy chain B, NMMHC-B
Immunogen	E. coli - derived recombinant Human MYH10 (His438-Pro716).
Accession	P35580
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.