

Product Information	
Catalog No.	PA10749
Product name	Anti-DDOST Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Dog, Rat, Bos taurus, Danio rerio, Chicken, Pongo abelii, Pig, Xenopus laevis, Xenopus tropicalis
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	DDOST 48 kDa subunit, OST48, Oligosaccharyl transferase 48 kDa subunit, KIAA0115, DDOST, Dolichyl-diphosphooligosaccharide--protein glycosyltransferase 48 kDa subunit
Immunogen	E. coli - derived recombinant Human DDOST (Ser43-Pro427).
Accession	P39656
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.