

Product Information	
Catalog No.	PA10752
Product name	Anti-CD130/IL6ST Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	IL-6 receptor subunit beta, CDw130, IL-6R-beta, Membrane glycoprotein 130, CD130, Oncostatin-M receptor subunit alpha, Interleukin-6 signal transducer, IL-6RB, IL-6R subunit beta, IL6ST, Interleukin-6 receptor subunit beta, gp130
Immunogen	E. coli - derived recombinant Human CD130/IL6ST (Thr392-Thr566).
Accession	P40189
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.