


Product Information	
Catalog No.	PA10758
Product name	Anti-STAT3 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat, Bos taurus, Chicken, Pig, Xenopus laevis
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Signal transducer and activator of transcription 3, APRF, STAT3, Acute-phase response factor
Immunogen	E. coli - derived recombinant Human STAT3 (Met1-Met770).
Accession	P40763
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-STAT3 polyclonal antibody at 1ug/mL on various samples.

Lane 1: Hela cell lysate

Lane 2: A549 cell lysate

Lane 3: A431 cell lysate

Lane 4: 293T cell lysate