

Product Information	
Catalog No.	PA10792
Product name	Anti-CHRNA4 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat, Chicken, Mustela putorius furo, Pan troglodytes
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Neuronal acetylcholine receptor subunit alpha-4, NACRA4, CHRNA4
Immunogen	E. coli - derived recombinant Human CHRNA4 (Asn49-Leu242).
Accession	P43681
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.