

Product Information	
Catalog No.	PA10821
Product name	Anti-FNTA Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat, Bos taurus
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Protein farnesyltransferase/geranylgeranyltransferase type-1 subunit alpha, Type I protein geranyl-geranyltransferase subunit alpha, FTase-alpha, FNTA, GGTase-I-alpha, CAAX farnesyltransferase subunit alpha, Ras proteins prenyltransferase subunit alpha
Immunogen	E. coli - derived recombinant Human FNTA (Met51-Lys366).
Accession	P49354
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.