

Product Information	
Catalog No.	PA10845
Product name	Anti-RORC/NR1F3/RORgT Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Pongo abelii
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	RAR-related orphan receptor C, RORC, NR1F3, Nuclear receptor ROR-gamma, Nuclear receptor RZR-gamma, Retinoid-related orphan receptor-gamma, RZRG, Nuclear receptor subfamily 1 group F member 3, RORG, RORgT
Immunogen	E. coli - derived recombinant Human RORC/NR1F3/RORgT (Glu212-Ser461).
Accession	P51449
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.