

Product Information	
Catalog No.	PA10918
Product name	Anti-RAB11A Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat, Bos taurus, Caenorhabditis elegans, Dog, Chlamydomonas reinhardtii, Diplobatis ommata, Chicken, Oryctolagus cuniculus, Pongo abelii, Pig
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Rab-11, RAB11A, RAB11, YL8, Ras-related protein Rab-11A
Immunogen	E. coli - derived recombinant Human RAB11A (Met1-Glu162).
Accession	P62491
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.