

Product Information	
Catalog No.	PA10967
Product name	Anti-TAGLN Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat, Bos taurus, Chicken
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	SM22-alpha, WS3-10, Smooth muscle protein 22-alpha, SM22, Transgelin, TAGLN, Protein WS3-10, 22 kDa actin-binding protein
Immunogen	E. coli - derived recombinant Human TAGLN (Met1-Ser201).
Accession	Q01995
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.