

Product Information	
Catalog No.	PA10973
Product name	Anti-GSTM4 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Mouse, Human, Macaca fascicularis, Rat
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	GSTM4-4, Glutathione S-transferase Mu 4, GSTM4, GST-Mu2, Leukotriene C4 synthase GSTM4, GST class-mu 4
Immunogen	E. coli - derived recombinant Human GSTM4 (Arg18-Trp215).
Accession	Q03013
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.