

| Product Information       |                                                                                                                                                                                 |
|---------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog No.               | PA11005                                                                                                                                                                         |
| Product name              | Anti-PPID Polyclonal Antibody                                                                                                                                                   |
| Applications              | ELISA, IHC, WB                                                                                                                                                                  |
| Species reactivity        | Bos taurus, Human, Mouse, Rat                                                                                                                                                   |
| Host species              | Rabbit                                                                                                                                                                          |
| Clonality                 | Polyclonal                                                                                                                                                                      |
| Isotype                   | IgG                                                                                                                                                                             |
| Form                      | Liquid                                                                                                                                                                          |
| Storage buffer            | 0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.                                                                                                                             |
| Purification              | Purified by antigen affinity column.                                                                                                                                            |
| Target                    | Rotamase D, 40 kDa peptidyl-prolyl cis-trans isomerase, PPID, Cyclophilin-40, Peptidyl-prolyl cis-trans isomerase D, PPIase D, CYP40, CYP-40, Cyclophilin-related protein, CYPD |
| Immunogen                 | E. coli - derived recombinant Human PPID (Gln7-Pro206).                                                                                                                         |
| Accession                 | Q08752                                                                                                                                                                          |
| Stability and Storage     | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.       |
| Product Usage Information | ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000                                                                                                                             |
| Note                      | For research use only.                                                                                                                                                          |