

Product Information	
Catalog No.	PA11081
Product name	Anti-CASP8 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	CAP4, MCH5, FADD-homologous ICE/ced-3-like protease, FLICE, Apoptotic cysteine protease, MORT1-associated ced-3 homolog, MACH, FADD-like ICE, ICE-like apoptotic protease 5, CASP8, CASP-8, Caspase-8, Apoptotic protease Mch-5
Immunogen	E. coli - derived recombinant Human CASP8 (Ser217-Asp479).
Accession	Q14790
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.