

Product Information	
Catalog No.	PA11156
Product name	Anti-SIGIRR Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	TIR8, Toll/interleukin-1 receptor 8, Single Ig IL-1R-related molecule, Single Ig IL-1-related receptor, Single immunoglobulin domain-containing IL1R-related protein, SIGIRR
Immunogen	E. coli - derived recombinant Human SIGIRR (Met1-His118).
Accession	Q6IA17
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.