


Product Information	
Catalog No.	PA11277
Product name	Anti-SEMA3C Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Bos taurus, Chicken, Pongo abelii
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	SEMA3C, Semaphorin-E, Semaphorin-3C, SEMAE, Sema E
Immunogen	E. coli - derived recombinant Human SEMA3C (Arg28-Asp150).
Accession	Q99985
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-SEMA3C polyclonal antibody at 1ug/ml on various samples.

- Lane 1: Raw264.7 cell lysate
- Lane 2: MCF7 cell lysate
- Lane 3: PC12 cell lysate
- Lane 4: Jurkat cell lysate
- Lane 5: Mouse lung lysate
- Lane 6: Mouse small intestine lysate
- Lane 7: Mouse heart lysate
- Lane 8: Rat lung lysate
- Lane 9: Rat small intestine lysate
- Lane 10: Rat heart lysate