

Product Information	
Catalog No.	PA11321
Product name	Anti-NANOG Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Macaca fascicularis, Pan paniscus, Pan troglodytes
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Homeobox protein NANOG, NANOG, Homeobox transcription factor Nanog, hNanog
Immunogen	E. coli - derived recombinant Human NANOG (Trp153-Val305).
Accession	Q9H9S0
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.