

Product Information	
Catalog No.	PA11343
Product name	Anti-ACP6 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Pongo abelii
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Acid phosphatase-like protein 1, ACPL1, PACPL1, ACP6, LPAP, Acid phosphatase 6, lysophosphatidic, Lysophosphatidic acid phosphatase type 6
Immunogen	E. coli - derived recombinant Human ACP6 (Met1-Glu428).
Accession	Q9NPH0
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.