

Product Information	
Catalog No.	PA11374
Product name	Anti-PDP1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat, Bos taurus, Pongo abelii
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	PDP 1, Pyruvate dehydrogenase phosphatase catalytic subunit 1, PDP1, PDPC 1, [Pyruvate dehydrogenase [acetyl-transferring]]-phosphatase 1, mitochondrial, PPM2C, Protein phosphatase 2C, PDP
Immunogen	E. coli - derived recombinant Human PDP1 (Asn315-Glu537).
Accession	Q9P0J1
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.