

Product Information	
Catalog No.	PA11404
Product name	Anti-APLN Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Homo sapiens, Bos taurus, Rattus norvegicus, Mus musculus
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	APEL, Apelin, APJ endogenous ligand, APLN
Immunogen	E. coli - derived recombinant Human APLN (Ser24-Phe77).
Accession	Q9ULZ1, Q9TUI9, Q9R0R3, Q9R0R4
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.