

Product Information	
Catalog No.	PA11458
Product name	Anti-EMC6/TMEM93 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Bos taurus, Danio rerio, Human, Mouse, Xenopus laevis, Xenopus tropicalis
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Transmembrane protein 93, TMEM93, EMC6, ER membrane protein complex subunit 6
Immunogen	E. coli - derived recombinant Human EMC6/TMEM93 (Val4-Thr47).
Accession	Q9BV81
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.