

Product Information	
Catalog No.	PA11466
Product name	Anti-Polyprotein pp62/Ba71V-94/ASFV-p15 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	African swine fever virus
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Polyprotein pp62, 60 kDa polyprotein, p60, 62 kDa polyprotein, p62, p15
Immunogen	E. coli - derived recombinant ASFV Polyprotein pp62/Ba71V-94/ASFV-p15 (Pro2-Asn160).
Accession	Q65179
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.