

Product Information	
Catalog No.	PA11476
Product name	Anti-YTHDC/LOC101245931 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Solanum chilense, Capsicum annuum, Solanum pennellii, Solanum tuberosum, Solanum lycopersicum, Anisodus tanguticus, Datura stramonium, Solanum dulcamara, Capsicum chinense, Solanum bulbocastanum, Solanum stoloniferum, Solanum pinnatisectum, Solanum stenotomum, Anisodus acutangulus, Solanum commersonii, Lycium ferocissimum, Lycium barbarum, Solanum verrucosum, Capsicum baccatum
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	YTH domain-containing protein, YTHDC, LOC101245931, Plant
Immunogen	E. coli - derived recombinant Tomato YTHDC/LOC101245931 (Met1-Ser215).
Accession	XP_004252938.1
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.