

Product Information	
Catalog No.	PA11486
Product name	Anti-NPPA Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Cat, Cavia porcellus
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Cardiodilatin, Pro atrial natriuretic factor 1-30, ANF 3-33, proANP, Atrial natriuretic factor prohormone, ANF 1-33, URO, LANH, proANP 31-67, Atrial natriuretic factor, Atrial natriuretic factor 8-33, VSDL, Atriopeptin I, CDD-ANF, ANF 8-33, Atrial natriuretic peptide prohormone, preproCDD-ANF, KP, CDD 95-126, Atrial natriuretic factor 3-33, Pro atrial natriuretic factor 99-126, Atrial natriuretic factor 1-33, proANP 95-126, proANF 99-126, proANF 79-98, Pro atrial natriuretic factor 31-67, ANP, Atriopeptin III, Alpha-hANP, Pro atrial natriuretic peptide 79-98, ANF, preproANP, Atriopeptin II, PND, Pro atrial natriuretic peptide 95-126, Natriuretic peptides A, Pro atrial natriuretic peptide 1-30, CDD-ANP (95-126), proANP 79-98, Pro atrial natriuretic factor 79-98, Atriopeptigen, proANF, CDD, proANF 1-30, LANP, CDD-ANP (99-126), Pro atrial natriuretic peptide 31-67, NPPA, Alpha-atrial natriuretic peptide, proANF 31-67, Cardionatrin, Long-acting natriuretic hormone, proANP 1-30
Immunogen	E. coli - derived recombinant Human NPPA (Ala25-Arg123).
Accession	P01160
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.