

Product Information	
Catalog No.	PA11515
Product name	Anti-CD158b2/KIR2DL3 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	NKAT2b, KIR2DL3, p58 natural killer cell receptor clone CL-6, p58.2 MHC class-I-specific NK receptor, CD158 antigen-like family member B2, MHC class I NK cell receptor, KIR-023GB, Killer inhibitory receptor cl 2-3, p58 NK receptor CL-6, CD158b2, Killer cell immunoglobulin-like receptor 2DL3, Natural killer-associated transcript 2, CD158B2, KIRCL23, NKAT-2, NKAT2a, NKAT2
Immunogen	E. coli - derived recombinant Human CD158b2/KIR2DL3 (His22-Pro341).
Accession	P43628
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.