

Product Information	
Catalog No.	PA11528
Product name	Anti-FKBP4 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat, Bos taurus, Oryctolagus cuniculus
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	PP1ase FKBP4, FKBP51, HSP-binding immunophilin, FKBP52, p59, HBI, Rotamase, Peptidyl-prolyl cis-trans isomerase FKBP4, FKBP4, FKBP59, 51 kDa FK506-binding protein, 52 kDa FK506-binding protein, FK506-binding protein 4, FKBP-4, 59 kDa immunophilin, FKBP-52, Immunophilin FKBP52, 52 kDa FKBP
Immunogen	E. coli - derived recombinant Human FKBP4 (Thr2-Ala459).
Accession	Q02790
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.