

Product Information	
Catalog No.	PA11554
Product name	Anti-NT5C3A Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Chicken
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Cytosolic 5'-nucleotidase III, P5N-1, 7-methylguanosine phosphate-specific 5'-nucleotidase, NT5C3, 7-methylguanosine nucleotidase, P5N1, p36, NT5C3A, cN-III, Cytosolic 5'-nucleotidase 3A, UMPH1, P5'N-1, Cytosolic 5'-nucleotidase 3, Uridine 5'-monophosphate hydrolase 1, Pyrimidine 5'-nucleotidase 1, PN-I
Immunogen	E. coli - derived recombinant Human NT5C3A (Met52-Leu336).
Accession	Q9H0P0
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.