

Product Information	
Catalog No.	PA11556
Product name	Anti-CD77/A4GALT Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Gorilla gorilla gorilla, Pan troglodytes, Pongo pygmaeus
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Gb3 synthase, UDP-galactose:beta-D-galactosyl-beta1-R 4-alpha-D-galactosyltransferase, CD77 synthase, Alpha-1, 4-N-acetylglucosaminyltransferase, Lactosylceramide 4-alpha-galactosyltransferase, Alpha4Gal-T1, Alpha-1, 4-galactosyltransferase, A4GALT1, A4GALT, Globotriaosylceramide synthase, P1/Pk synthase, A14GALT
Immunogen	E. coli - derived recombinant Human CD77/A4GALT (Gly44-Leu353).
Accession	Q9NPC4
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.