

Product Information	
Catalog No.	PA11567
Product name	Anti-SLEV NS1/Non-structural protein 1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	St. louis encephalitis virus
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	polyprotein, POLY, NS1, Non-structural protein 1
Immunogen	E. coli - derived recombinant St. Louis encephalitis virus (SLEV) NS1/Non-structural protein 1 (Asp790-Ala1141).
Accession	YP_001008348.1
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 4°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.