

Product Information	
Catalog No.	PA11607
Product name	Anti-Monkeypox virus/MPXV D13L Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Monkeypox virus, Vaccinia virus, Variola virus
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	D13L
Immunogen	E. coli - derived recombinant Monkeypox virus/MPXV D13L (Met1-Glu315).
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.